



THUASNE


Jump in the
softness



VENOFLEX

KOKOON ABSOLU

Medical compression
socks, stockings or tights
for women



1 in 2 women suffer from venous insufficiency*

They are looking for **premium solution** to meet 3 needs:**


QUALITY
Feel safe


AESTHETIC
Feel pretty


COMFORT
Feel at ease

A softness and casual phlebology range approved by patients with **visible premium attributes to improve compliance!****


COMFORT & INNOVATIVE FINISHES



A new simple fold-down on socks for a better holding without tourniquet effect



A new breathable antislip thanks to wavy silicone, to stay in place all day long



An extensible panty on tights without pressure on stomach

SOFTNESS

- Instantly & highly recognized by all users**
- To take care of patient legs, by the careful choice of yarns and specific construction

“ I would feel (...) comfortable in it. I would also have a sense of being taken care of. Someone is taking care of me ”

AESTHETIC

- Modern finishes: with Venoflex drop signature on all products
- For a casual style

“ We are moving away from the medical and getting closer to femininity and fashion ”

QUALITY

- Fully perceived**
- Stitching details
- Reinforced foot areas

“ Good quality material (...) because it does not catch, doesn't stretch permanently, they go back to their initial form ”

Thuasne gives **women** back their **confidence**, thanks to **KOKOON ABSOLU** they feel feminine, and **together** we **take care** of them.**

Did you know?

Compliance depends on the quality of information provided by pharmacists***


To support you, Thuasne offers a complete range of products to meet most patients' needs, whatever the season, morphology or clothing style. Use the decision tree to help you deliver the product. You are one of the core partners in compliance, help women to become actors in their own health.

*HAS**Focus group led in France + Poland - 29 patients - OPINION WAY - Oct 20 *** Gillet JL et al. Study in pharmacies of the determinants of the adherence to a medical prescription for compression stockings and of the people satisfaction. Phlébologie 2013, 66, 1, p.14-21

VENOFLEX

KOKOON ABSOLU

A wide range available in class 2 (15-20 mmHg) to take care of women's legs in all femininity.



Socks				
	cB (cm)	cC - (cm)	cC (cm)	cC + (cm)
0	16 - 19	21 - 28	28 - 33	33 - 42
1	19 - 21	23 - 30	30 - 35	35 - 46
2	21 - 24	26 - 32	32 - 38	38 - 48
3	24 - 26	29 - 34	34 - 42	42 - 52
4	26 - 29	32 - 37	37 - 46	46 - 54
5	29 - 31	36 - 43	43 - 49	49 - 56

Normal $dD < 40$ cm - Long $dD > 40$ cm

Stocking and tights				
	cB (cm)	cG - (cm) STOCKINGS ONLY	cG (cm)	cG + (cm) STOCKINGS ONLY
0	16 - 19	35 - 42	42 - 52	52 - 57
1	19 - 21	38 - 45	45 - 54	54 - 59
2	21 - 24	43 - 50	50 - 59	59 - 66
3	24 - 26	48 - 56	56 - 66	66 - 70
4	26 - 29	51 - 58	58 - 69	69 - 74
5	29 - 31	55 - 63	63 - 74	74 - 78

Normal $dK < 74$ cm - Long $dK > 74$ cm

Indications: 15-20 mmHg/Class 2: heavy legs, restless legs, evening oedema, telangiectasia or reticular veins (< 3 mm), following sclerotherapy for telangiectasia or reticular veins (< 3 mm), following sclerotherapy or surgery for varicose veins (≥ 3 mm), varicose veins (≥ 3 mm), prevention of venous thrombosis in pregnant women or post-partum, prevention of venous thrombosis during long journeys, treatment of lymphoedema (maintenance phase).

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region. All the products mentioned on this document are class I medical devices, CE marked according to the European Regulation 2017/745 on medical devices. Please read carefully the instructions for use, indications and contraindications of the product. 2104159 (2021-06) - Photo credit: Studio Caterin