



Automatic Manosonic Nebulizing System

Treatment of tubo-tympanic pathologies:

- secretory otitis media or seromucous otitis,
- tubotympanic sequelae (like retraction pocket...).

3 combined functions, automatically synchronized with the patient's swallowing, allow the aerosol penetration into the Eustachian Tube:

- aerosol production,
- emission of 100 Hz sonic vibrations,
- temporary application of a modulable adjustable overpressure.


Contra-indications:

- child under 4 years old,
- patient with a patulous Eustachian tube,
- patient with cleft lip and palate.




Latex free accessories: **PURENEB® EXTRANEB MANOSONIC** nebulizing kits with overpressure tight nose pieces.



 **EXTRANEB MANOSONIC** adult, ref. MS1A
medium and large nose pieces



 **EXTRANEB MANOSONIC** paediatric, ref. MS1E
small and medium nose pieces



Carrying bag, ref. F64



ATOMISOR CONTROL¹ to check the good functioning of your device, ref. MEA1

Class IIa Medical Devices. Carefully read the instructions for use provided with the medical device. Documentation for Health Professionals. Non contractual pictures



A. COMPS.EN.044.0120

AEROSOL THERAPY BY DTF MEDICAL



La Diffusion Technique Française
19, rue de la Presse - CS 60132 - 42003 Saint-Etienne Cedex 1 - France

Tél. : (33) 4 77 74 51 11 - Fax : (33) 4 77 79 67 72

E-mail : dtf@dtf.fr - www.dtf.fr

