

BENOIT SYSTEMES



MOBILI KIT



WWW.BENOITSYSTEMES.COM



MOBILI KIT

“Your seat becomes mobile !”

Do you want to preserve your independence ?

Do you want to be more active and efficient in everyday life?

MOBILI KIT is made for you.

It will help you to move around, both at home and at work.

Compact, simple and silent, this device knows how to become discreet under a seat.



SILVERSHOW AWARD

1ST PRIZE FOR THE BEST SOLUTION OF AGEING WELL AND LOSS OF AUTONOMY

CONCOURS LÉPINE AWARD


1ST PRIZE

HOW DOES IT WORK ?

1

2

3



CHOOSE YOUR FAVORITE
CHAIR


SET UP BY YOURSELF OR BY
A TECHNICIAN

TAKE A SEAT !

4

5

6



LIFT UP THE FOOTRESTS

MAKE YOUR LIFE EASIER!

ENJOY IT !



MOBILI KIT can be installed on any
type of 4-legged seat,
even your favourite armchair!

The joystick control system offers
great precision and manoeuvrability.



MOBILI KIT gets
started fast and easily.



“ **MOBILI KIT,**
Let's set off a
better active life ! ♡♡

MOBILI KIT ADVANTAGES



ADAPTS TO 80%
OF THE SEAT
WITH 4 LEGS



INSTALLATION BY
THE USER OR
BY A TECHNICIAN



SIMPLE AND
ACCURATE
OPERATED



MASTER IT IN
30 MIN



ELECTRIC TILT
SYSTEM



5 ADJUSTABLE
SPEED




4 DAYS
AUTONOMY



COMPACT

Its **compact design** allows you to move around, even into the tightest areas.



MOBILI KIT
5 times more active
with 5 times less efforts

THEY ARE CONVINCED, SO WHY NOT YOU ?

SYLVIE

User since 2019

Her mother Danielle,
tells us :

“ Sylvie used a walker to move in our house; but it was not safe enough. Today **Mobili Kit** avoids any fall: I feel more serene, and Sylvie enjoys full freedom of movement. She can go everywhere in the house and even go out on the terrace! For indoors, it is more suitable than a wheelchair because it is less bulky, especially for the corridors. We are very happy with it!



EMMANUEL

User since 2021

His wife Véronique,
tells us :

“ My husband was using a walker, but he had more and more difficulties, so we decided to purchase a **Mobili Kit** after meeting Robert Benoit. We are using it every day, anywhere in the house because the chair can spin around and pass the doors without having to modify their width. It is very convenient and comforting and it prevents you from falling!



Find more testimonials on the website
benoitsystemes.com

TECHNICAL DETAILS

Min. dimensions	L58 x W41 x H23 cm
Max. weight of user	130 kg
MOBILI KIT weight	20 kg
Adjustable speed	from 1 to 4 km/h (possibility to limit the speed by programming)
Max. safe slope	3° (around 5%)
Max. crossable obstacle height	2 cm


Possible dimensions of the platform
based on chosen chair

from 25 to 75 cm at the front,
from 25 to 68 cm at the rear,
from 32 to 67 cm in deep

Guarantee — 24 months for the motor unit and battery

BATTERY

Lithium-ion battery	9Ah - 227 Wh
Weight	3 kg
Autonomy	4 days





Conception graphique © E-to LINK - Photos © Philippe Prof / © Lionel Dupouy - 11/2022

REQUEST FOR A FREE TRIAL !

📍 7 Rue du Pont,
21450 Billy-Lès-Chanceaux
☎️ FRANCE (+33) 3 80 96 51 25
@ export@benoisystemes.com

CE MOBILI KIT complies with the European Directive 93/42/CEE on medical devices.